



Stark Genotyping HPV
Molecular Diagnostic Kit



Contents

Description.....	3
Intended Use	3
Kit Contents.....	3
Storage	3
Guarantee and Warranty	3
Warning and Precautions	3
Quality Control	4
Materials Required (but Not Provided)	4
Procedures	4
Applications	5
Features.....	5
Recommended Starting Material	5
Protocol	7
Results Interpretation	7
Control Conditions for a Valid PCR Run	8
Test Limitations.....	9
Performance Evaluation.....	9
Symbols	40
Troubleshooting	40
Technical assistance	41
Contact Information	41

Description

The Stark Genotyping HPV Molecular Diagnostic Kit utilizes in vitro nucleic acid TaqMan assay technology, employing polymerase chain reaction for the genotyping detection of 14 high-risk and two low-risk types of HPV DNA in cervical or vaginal specimens. The high-risk HPV types detected include 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68, alongside low-risk types 6 and 11. This test is intended exclusively for professional use by trained and validated laboratory personnel. The Kit employs a three-Channels Multiplex Real-Time PCR system and includes an internal control to enhance isolation quality and prevent false PCR inhibition results.

Intended Use

The Stark HPV Genotyping Molecular Diagnostic Kit (Real-Time PCR) is intended for in vitro diagnostic use. This Kit is designed for the qualitative detection and genotyping of human papillomavirus (HPV) DNA in cervical or vaginal specimens. Results obtained are for professional use only by trained and validated laboratory personnel. This Kit aids in identifying HPV infections, assisting healthcare professionals in making informed decisions regarding patient management and treatment strategies. Careful adherence to provided instructions is crucial to ensure accurate and reliable results.

Kit Contents

Ingredients	25 Preps (REF: ST242001)
Q-ROMAX, 4X	600µl
Pro 1 HPV Mix	350µl
Pro 2 HPV Mix	350µl
Pro 3 HPV Mix	350µl
Pro 4 HPV Mix	350µl
Positive Control	150µl
Water (PCR Grade)	150µl

Storage

Ensure all reagents are stored away from direct light in darkness at temperatures ranging between -20°C to -25°C. Avoid frequent freeze-thaw cycles. When stored as directed, all reagents maintain stability until the expiration date indicated on the Kit box.

Guarantee and Warranty

CARBON Technologies LLC assures the efficiency of all manufactured Kits and reagents. For guidance on selecting the appropriate Kits for your requirements, please reach out to our technical support team. Should the products fail to meet your satisfaction for reasons other than misuse, please contact our technical support team. In the event of issues stemming from the manufacturing process, CARBON Technologies LLC will promptly replace the Kit.

Warning and Precautions

- This Kit is for in vitro diagnostic use only.
- Material Safety Data Sheets (MSDS) for all products and reagents are available online at www.carbontechnologiesco.com.
- Adhere to laboratory safety protocols diligently.
- Familiarize yourself thoroughly with the guidelines before usage.
- Prohibit eating, drinking, smoking, chewing gum, applying cosmetics, or taking medicine in laboratories where hazardous materials and human samples are handled.
- Treat all patient samples and positive controls as potentially infectious.

- Use the Stark Genotyping HPV Molecular Diagnostic Kit under the supervision of a physician for emergency and in vitro diagnostic purposes.
- Ensure all steps of procedures, including sampling, storage, shipping, and laboratory tests, comply with biosafety and molecular laboratory management standards.
- Equip clinical laboratories with instruments and operators in accordance with the regulations of the Ministry of Health.
- Any alteration or replacement of Kit contents may affect functionality and contravene product licensing.
- Utilize sterile and DNase-RNase-free pipette tips and microtubes to prevent contamination. Change filter pipette tips after each substance or sample addition.
- Dispose of waste in line with biosafety guidelines. Regularly sanitize desks and laboratory instruments with 70% Ethanol or 10% Sodium Hypochlorite.
- Shield Pro Mixes from sunlight exposure.
- Clean and disinfect specimen spills promptly using suitable disinfectants adhering to national and local regulations.
- Dispose of all specimens, reagents, and potentially contaminated materials following national and local regulations.
- The hazard and precautionary statements provided pertain to the components of the Stark Genotyping HPV Molecular Diagnostic Kit.

Quality Control

The Stark Genotyping HPV Molecular Diagnostic Kit undergoes rigorous testing against predetermined experiments on a lot-to-lot basis to ensure consistent product quality. Accessible results of these experiments are obtainable online by referencing the REF and Lot numbers at www.carbontechnologiesco.com.

Materials Required (but Not Provided)

- DNase-RNase-free microtubes (1.5ml)
- PCR microtube 0.1ml or 0.2ml strip
- Various models of pipettes and pipette tips (10 μ l, 100 μ l, and 1000 μ l filter pipette tips)
- Surface sanitizing solution
- Disposable Powder-Free gloves and surgical gown
- Three-Channels Multiplex Real-Time PCR Instrument (with green, yellow, and orange channels)
- Vortex
- Cool box

Procedures

The Stark Genotyping HPV Molecular Diagnostic Kit employs a polymerase chain reaction of Real-Time PCR for its testing procedure. Specifically designed for genotyping molecular diagnosis of L1, E1, E2, E6, and E7 genes of the Human Papillomavirus (HPV), this Kit facilitates nucleic acid isolation using the DNall VirAll Kit or other approved Kits. Verified sample combinations are then added to the master mix primer/probe mix to initiate the reaction. Additionally, the Stark Genotyping HPV Molecular Diagnostic Kit incorporates a second heterologous amplification system to identify potential PCR inhibition. This is detected as an internal control (IC) in the fluorescence channel Cycling Yellow of Real-Time PCR instruments. Through its sampling mechanism, the quality of sample isolation and PCR reaction process can be monitored and controlled to prevent false-negative results. The assay demonstrates a Limit of Detection (LoD) of 5 Copies/5 μ l for types 16, 18, and 45, 50 Copies/5 μ l for types 6/11, 51, 56/66, 33/52/58, 35, 59, and 39, and 500 Copies/5 μ l for types 68 and 31.

Applications

The Stark Genotyping HPV Molecular Diagnostic Kit technology is an in vitro nucleic acid TaqMan assay employing signal amplification through polymerase chain reaction and fluorescent probes (ROX/Texas Red, Yakima Yellow, and FAM). It is utilized for the genotyping detection of 14 high-risk and two low-risk types of HPV DNA in cervical or vaginal specimens. The high-risk HPV types detected include 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68, along with low-risk types 6 and 11. This diagnostic test Kit is applicable in human cervical screening, cytology samples, urine, and paraffin-embedded tissue analysis. Utilizing the polymerase chain reaction (PCR), the Kit is configured with three-Channels Multiplex Real-Time PCR instruments.

Features

Table 1: Stark Genotyping HPV Molecular Diagnostic Kit features and specifications

Technology	Real-Time PCR
Type of Analysis	Genotyping
Target Sequence	L1, E1, E2, E6 and E7 genes
Analytical Specificity	high-risk HPV types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 low-risk HPV types 6 and 11
Analytical Sensitivity (LOD)	for types 16, 18, 45 is 5 Copies/5µl, for types 6/11, 51, 56/66, 33/52/58, 35, 59, 39 is 50 Copies/5µl and for types 68, 31 is 500 Copies/5µl.
Diagnostic Specificity	100%
Diagnostic Sensitivity	98.52%
Extraction/Inhibition Control	PCR inhibition and DNA extraction efficiency control
Validated Specimen	Cervical screening, Pap smear, Urine and Paraffin-embedded tissue sample
Storage	-20 ± 5°C
Validated Extraction Method	DNall VirAll Kit
Instruments	Rotor-Gene Q, 2plex, Corbett Rotor-Gene 3000&6000, Mic qPCR Cycler, StepOne and StepOne plus Applied Biosystem
Required Detection Channels	Green-Yellow-Orange

Recommended Starting Material

- Before commencing any tests, ensure each component is thawed, vortexed, and briefly centrifuged. Avoid subjecting the components to repeated freeze-thaw cycles.
- Treat all samples as infectious and biohazardous, following safe laboratory procedures.
- Collect samples from the cervix using a brush or swab, cytology samples, urine, or paraffin-embedded tissue.

Storage and Sample Preparation

Cervical screening, Pap smear, Urine, and Paraffin-embedded tissue samples:

Fresh specimens should either be processed immediately according to the sample procedure outlined in the Sample Processing Protocol section or stored frozen at -20°C. Frozen samples must be brought to room temperature before initiating sample processing. Sample Pre-treatment decontaminates the specimen and prepares it for extraction.

DNA Isolation

For nucleic acid isolation, utilize the DNall VirAll Kit or other approved Kits.

Before Starting

Remove each component from the Kit and place them on the benchtop. Allow the reagents to equilibrate to room temperature, then briefly vortex each tube for later use.

Buffer Preparation

Take out each component from the Kit and allow the reagents to equilibrate to room temperature. Before use, briefly vortex the components. The total volume of isolated nucleic acid should be 5µl. Refer to Table 1 for buffer preparation and Table 2 for PCR run instructions.

Table 2: Preparation of components per single reaction

Components	Volume
Q-ROMAX, 4X	6µl
Pro 1 HPV or Pro 2 HPV or Pro 3 HPV or Pro 4 HPV Mix	14µl
Isolated DNA	5µl

Pathogenicity

Human papillomaviruses (HPVs) are small double-stranded DNA viruses that infect the cutaneous and mucosal epithelium. Infection by specific HPV types has been associated with the development of cervical carcinoma. HPV targets epithelial cells undergoing terminal differentiation and employs multiple mechanisms to override normal differentiation regulation, producing progeny virions. This leads to morphological changes resulting in noncancerous tumors in the skin, commonly referred to as "Papilloma".

Human papillomavirus comprises 150 identified strains categorized into various types. Approximately 75% of papilloma types target epithelial cells (cutaneous types), while the remaining 25% (around four types) target mucosal epithelial tissues, identified as mucosal or genital types. Cutaneous types mainly cause noncancerous epithelial warts in different skin areas, particularly in the lips' epithelial cells, while mucosal types induce genital warts and cancerous tumors, especially cervical cancer.

Low-risk HPV

Low-risk HPV strains include types that cause noncancerous and low-risk warts in genital areas of both women and men. In women, they predominantly manifest in internal genitalia, such as the vagina and cervix. There are 12 types of papillomaviruses with low cancer risk that can induce warts and morphological changes in the genital area: 6, 11, 40, 42, 43, 44, 53, 54, 61, 72, 73, and 81. Types 6 and 11 are responsible for 90% of genital warts.

High-risk HPV

High-risk HPV can lead to cancers of the cervix, vagina, and vulva in women; penis in men; anus; and the back of the throat, including the base of the tongue and tonsils (oropharynx), in both men and women. There are 14 types of high-risk HPV viruses, including: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68. Approximately 70% of reported cervix cancers are caused by HPV 16 and 18. Additionally, HPV 16 is responsible for 95% of oral, vaginal, and anal cancers. HPV 31 and 45 are considered high-risk cancer strains, contributing to 5%-10% of cervix cancers.

Workstation Preparation

Before usage, all work surfaces, pipettes, and other supplies must be thoroughly cleaned and sanitized. To minimize

the risk of nucleic acid contamination, utilize sanitizers such as 70% Ethanol or 10% Sodium Hypochlorite.

Protocol




Figure 1: preparation of reagents, PCR run, and interpretation of results.

PCR reaction preparation

Table 3: PCR reaction preparation

Components	Volume
Q-ROMAX, 4X	6µl
Pro1 HPV or Pro2 HPV or Pro3 HPV or Pro4 HPV or Mix	14µl
Isolated DNA	5µl

Thermal Profile

Table 4: Thermal profile for Stark Genotyping HPV Molecular Diagnostic Kit

Cycles	Temperature	Incubation Time	Cycle Numbers
Pre-Denaturation	95 °C	3 min	1
Denaturation	95 °C	10 sec	45
Annealing and acquisition on channel Green, Yellow and Orange	57°C	30 sec	

Results Interpretation

Data analysis for each type should be performed separately using a manual threshold. For the detection of L1, E1, E2, E6, and E7 genes, fluorophores FAM (green), Yakima Yellow (Yellow), and ROX/Texas Red (orange) are utilized for Pro1 HPV Mix to Pro4 HPV Mix. The IC gene in Pro2 HPV Mix is detected using the Yakima Yellow (Yellow) fluorophore.

Table 5: Specific fluorescent channels identified for HPV types.

ProMIX	Green	Yellow	Orange
Pro1 HPV	45	16	18
Pro2 HPV	51	Internal Control	56/66
Pro3 HPV	35/39	33/52/58	6/11
Pro4 HPV	31	68	59

Control Conditions for a Valid PCR Run

A negative control serves as a contamination control. If the magnitude increase of the fluorescence curve in the negative control does not cross the threshold and Ct is less than 35 (Ct<35), it is considered as possible contamination. Strong signals above 35 in the NTC can be PCR artifacts. In such cases, the shape of the curve, typically an S-shaped curve, should be considered for a positive result.

The internal control should yield a positive result for all clinical specimens with Ct ≤ 35, indicating sufficient nucleic acid from the human gene and acceptable sample quality. If the internal control curve has Ct > 37 or lacks Ct, it indicates low sample concentration or inhibitors in the reaction. In such instances, diluting the isolated sample by at least half is recommended. If the test result remains unacceptable upon retesting, obtain a new sample from the patient and repeat the test.

A positive clinical specimen should have Ct ≤ 37 for the gene. If the expected positive reaction, characterized by a typical S-shaped curve, is not achieved, the performed test is deemed unacceptable. Repeat the test following the instructions provided in the Kit catalog.

In case of failure of the positive control, determine the reason, take corrective action, and document the results of the corrective action.

Table 6: Control conditions for a valid PCR Run

ProMIX	ROX/Texas Red	Yakima Yellow	FAM	Results
Pro 1 HPV	+	+	+	Positive:45 Positive: 16 Positive:18
Pro 2 HPV	+	It is not considered	+	Positive:51 Positive:56/65
Pro 3 HPV	+	+	+	Positive:35/39 Positive:33/52/58 Positive:6/11
Pro 4 HPV	+	+	+	Positive:31 Positive:68 Positive:59
Pro 2 HPV	-	+	-	Negative Result
Pro 1 HPV or Pro 2 HPV or Pro 3 HPV or Pro 4 HPV	-	-	-	Invalid results

Test Limitations

- Low virus titers in patients' specimens, improper transportation, and low-quality DNA isolation can lead to false-negative results.
- All related controls should be checked before result interpretation to ensure reliable results. Failure to do so may compromise the accuracy of the results.
- The limit of detection (LoD) of the present Kit is demonstrated by Ct ≤ 37, and a typical S-shaped curve must appear for all positive specimens.
- Improper storage conditions can result in false-negative results.
- The product should only be used by personnel specially instructed and trained in in-vitro diagnostics, as individual errors can compromise the accuracy of the results.
- Diagnosis of Human papillomavirus infection is made when test results are positive and accompanied by clinical symptoms. Treatment should be conducted based on diagnostic Kit results, medical background, and response to treatment.

Performance Evaluation

- **Limit of Detection (LoD) - Analytical Sensitivity**

LoD studies were conducted to determine the lowest detectable concentration of HPV DNA, at which approximately 95% of all true positive replicates test positive. The LoD was determined through limiting dilution studies using characterized samples. Analytical sensitivity, considering the purification (DNall VirAll Kit) of the Stark Genotyping HPV Molecular Diagnostic Kit, was determined using a dilution series of standards ranging from 5 to nominal 50 and 500 HPV genome equivalents (GE)/5µl spiked in clinical cervical specimens.

The LoD of each test was confirmed by testing 20 replicates with a dilution series (500, 50, 5 genome equivalents (GE)/5µl) at the tentative limit of detection. The final LoD of each test was determined as the lowest dilution series resulting in positive detection of 19 out of 20 replicates.

The LoD of the Stark Genotyping HPV Molecular Diagnostic Kit was established using the DNall VirAll Kit. Results demonstrated that the LoD of the assay for types 16, 18, and 45 is 5 Copies/5µl, for types 6/11, 51, 56/66, 33/52/58, 35, 59, and 39 is 50 Copies/5µl, and for types 68 and 31 is 500 Copies/5µl.

Table7: Detection Results of Stark Genotyping HPV Molecular Diagnostic Kit Using DNall VirAll Kit.

Test No	(Copies/5µl)		
	HPV type 16		
	500	50	5
1	32.14	34.04	36.04
2	32.51	35.14	36.14
3	32.57	35.37	36.47
4	32.01	34.57	35.41
5	32.11	34.91	34.95
6	31.75	34.87	35.21
7	32.67	34.18	36.74
8	32.84	35.77	37.54
9	31.47	35.24	36.25

10	32.55	34.44	37.64
11	31.94	34.15	Undetermined
12	31.27	34.27	35.73
13	31.53	35.04	36.28
14	31.84	35.21	37.42
15	31.95	35.11	36.08
16	32.01	35.29	36.47
17	32.43	34.77	37.51
18	31.47	34.69	37.69
19	32.63	34.94	36.73
20	32.77	34.23	36.47
Positive percentage in each concentration	100%	100%	95%

Test No	(Copies/5µl)		
	HPV type 18		
	500	50	5
1	33.19	35.18	37.10
2	34.33	35.46	36.87
3	33.61	36.35	36.49
4	33.12	36.49	Undetermined
5	34.40	35.13	36.77
6	33.07	35.02	37.14
7	33.41	35.49	37.04
8	34.34	36.20	37.24
9	33.53	35.82	38.08
10	33.89	35.61	37.53
11	34.06	36.63	36.94
12	33.38	36.36	38.91
13	33.29	35.98	36.28
14	34.72	36.62	37.10
15	33.76	35.15	37.27
16	34.79	35.11	37.24
17	33.78	36.46	37.91
18	33.22	35.53	37.21
19	34.44	36.06	37.66
20	33.69	36.15	37.14
Positive percentage in each concentration	100%	100%	95%

Test No	(Copies/5µl)		
	HPV type 45		
	500	50	5

1	33.19	35.18	37.10
2	34.33	35.46	36.87
3	33.61	36.35	36.49
4	33.12	36.49	Undetermined
5	34.40	35.13	36.77
6	33.07	35.02	37.14
7	33.41	35.49	37.04
8	34.34	36.20	37.24
9	33.53	35.82	38.08
10	33.89	35.61	37.53
11	34.06	36.63	36.94
12	33.38	36.36	38.91
13	33.29	35.98	36.28
14	34.72	36.62	37.10
15	33.76	35.15	37.27
16	34.79	35.11	37.24
17	33.78	36.46	37.91
18	33.22	35.53	37.21
19	34.44	36.06	37.66
20	33.69	36.15	37.14
Positive percentage in each concentration	100%	100%	95%

Test No	(Copies/5µl)		
	HPV type 33/52/58		
	500	50	5
1	36.17	37.16	Undetermined
2	35.42	38.06	39.38
3	35.61	37.61	39.49
4	35.52	37.82	39.52
5	35.49	38.13	41.01
6	36.07	38.41	41.48
7	35.41	37.49	41.59
8	35.37	37.84	Undetermined
9	35.53	37.71	Undetermined
10	35.89	37.61	Undetermined
11	36.06	37.88	Undetermined
12	35.36	37.46	40.91
13	35.52	38.02	41.28
14	35.92	37.49	41.60
15	35.76	37.92	Undetermined
16	35.62	37.91	Undetermined
17	36.71	37.77	Undetermined

Stark Genotyping HPV Molecular Diagnostic Kit

CARBONTECHNOLOGIES

18	36.47	37.15	Undetermined
19	36.51	38.08	Undetermined
20	36.01	38.25	Undetermined
Positive percentage in each concentration	100%	100%	45%

Test No	(Copies/5µl)		
	HPV type 39		
	500	50	5
1	36.61	38.11	Undetermined
2	37.37	Undetermined	Undetermined
3	37.88	38.15	Undetermined
4	37.02	38.19	Undetermined
5	37.59	38.72	Undetermined
6	37.07	38.02	Undetermined
7	37.41	37.82	41.04
8	37.92	37.60	Undetermined
9	36.53	37.42	41.58
10	36.89	37.72	40.83
11	37.06	37.53	Undetermined
12	37.41	37.73	40.91
13	37.28	37.92	41.48
14	37.01	38.16	Undetermined
15	36.03	38.35	Undetermined
16	37.38	37.91	Undetermined
17	36.79	38.12	Undetermined
18	37.12	38.02	41.31
19	36.66	38.27	Undetermined
20	37.50	37.85	Undetermined
Positive percentage in each concentration	100%	95%	30%

Test No	(Copies/5µl)		
	HPV type 59		
	500	50	5
1	Undetermined	37.78	40.10
2	35.44	38.06	40.17
3	35.83	37.66	41.59
4	36.02	37.92	40.58
5	35.89	Undetermined	40.44
6	36.07	37.02	Undetermined
7	35.48	37.19	41.04
8	35.55	38.60	Undetermined

9	35.53	38.12	40.58
10	35.89	37.51	41.53
11	36.06	38.12	Undetermined
12	36.16	38.66	41.91
13	36.28	37.98	42.28
14	35.92	37.66	40.60
15	35.16	38.15	Undetermined
16	35.76	39.94	Undetermined
17	35.66	38.46	Undetermined
18	36.22	39.09	Undetermined
19	36.04	37.88	40.15
20	36.44	37.52	Undetermined
Positive percentage in each concentration	95%	95%	60%

Test No	(Copies/5µl)		
	HPV type 35		
	500	50	5
1	36.27	38.28	41.01
2	37.20	39.06	41.52
3	37.61	38.15	41.49
4	36.02	39.79	37.04
5	Undetermined	39.13	39.22
6	37.07	40.02	Undetermined
7	37.40	40.49	39.04
8	37.37	40.60	Undetermined
9	36.53	39.42	40.58
10	37.89	Undetermined	Undetermined
11	38.06	40.63	Undetermined
12	37.36	39.86	41.21
13	38.18	39.98	Undetermined
14	36.92	39.56	40.10
15	37.76	40.15	Undetermined
16	36.66	41.01	Undetermined
17	37.72	40.46	Undetermined
18	36.21	39.33	Undetermined
19	37.44	39.06	Undetermined
20	37.31	39.09	Undetermined
Positive percentage in each concentration	95%	95%	45%

Test No	(Copies/5µl)		
	HPV type 68		

	500	50	5
1	37.27	Undetermined	Undetermined
2	37.12	Undetermined	Undetermined
3	37.71	40.13	Undetermined
4	38.92	40.49	Undetermined
5	Undetermined	39.81	Undetermined
6	39.07	40.21	Undetermined
7	37.41	Undetermined	Undetermined
8	38.37	39.21	Undetermined
9	37.13	40.12	Undetermined
10	38.89	Undetermined	Undetermined
11	38.12	Undetermined	Undetermined
12	37.31	39.56	40.91
13	38.28	40.18	40.32
14	38.92	Undetermined	41.04
15	38.16	Undetermined	Undetermined
16	37.23	Undetermined	Undetermined
17	37.72	Undetermined	Undetermined
18	38.12	Undetermined	39.78
19	38.54	Undetermined	39.17
20	38.29	Undetermined	Undetermined
Positive percentage in each concentration	95%	40%	2%5

Test No	(Copies/5µl)		
	HPV type 6/11		
	500	50	5
1	35.27	37.28	Undetermined
2	35.87	38.06	41.57
3	35.61	37.84	41.49
4	36.02	37.14	Undetermined
5	35.49	Undetermined	Undetermined
6	37.07	37.02	Undetermined
7	35.41	37.49	40.04
8	35.37	38.60	40.24
9	36.53	37.42	40.58
10	35.89	38.47	41.47
11	37.06	39.62	Undetermined
12	35.16	36.46	40.41
13	35.28	38.48	41.28
14	35.92	38.66	40.60
15	36.76	38.17	40.17
16	36.12	37.91	Undetermined

17	36.04	38.74	41.28
18	37.01	37.19	41.41
19	36.57	37.11	40.65
20	36.77	37.25	Undetermined
Positive percentage in each concentration	100%	95%	65

Test No	(Copies/5µl)		
	HPV type 66/56		
	500	50	5
1	39.71	Undetermined	Undetermined
2	39.30	39.06	Undetermined
3	39.41	39.35	Undetermined
4	37.02	40.49	Undetermined
5	38.49	41.13	41.34
6	38.64	41.02	Undetermined
7	37.14	40.49	42.05
8	37.37	39.60	Undetermined
9	38.53	39.42	40.28
10	37.89	40.61	Undetermined
11	38.06	39.48	Undetermined
12	37.36	41.46	40.91
13	37.78	40.48	40.28
14	33.92	40.66	40.60
15	34.76	39.15	Undetermined
16	34.76	39.97	40.38
17	35.72	39.47	Undetermined
18	38.71	39.41	40.25
19	37.97	38.67	41.25
20	37.57	39.24	Undetermined
Positive percentage in each concentration	100%	95%	45%

Test No	(Copies/5µl)		
	HPV type 51		
	500	50	5
1	34.73	39.28	Undetermined
2	34.57	39.06	Undetermined
3	35.21	40.35	Undetermined
4	35.47	39.27	Undetermined
5	34.89	39.13	Undetermined
6	36.07	39.27	Undetermined
7	35.41	39.57	41.37

8	35.37	39.97	40.97
9	34.98	Undetermined	40.58
10	35.74	39.61	40.53
11	35.57	39.62	Undetermined
12	36.36	39.67	40.91
13	35.28	39.17	Undetermined
14	34.92	39.69	39.60
15	34.73	39.18	Undetermined
16	34.47	39.91	Undetermined
17	35.01	40.46	39.27
18	37.27	40.19	40.41
19	36.12	39.57	40.65
20	35.74	39.95	Undetermined
Positive percentage in each concentration	100%	95%	45%

- Inclusivity (Analytical Sensitivity)**

The inclusivity of the primer/probe set utilized in the Stark Genotyping HPV Molecular Diagnostic Kit was assessed in silico using HPV sequences obtained from the NCBI database, accessed on September 26, 2021. Sequence alignment analysis of the primer/probe sets targeting the L1, E1, E2, E6, and E7 genes demonstrated 100% inclusivity for HPV sequences identified from patient samples. The table below illustrates representative alignment results for L1, E1, E2, E6, and E7 genes.

Table 8: Alignment test result for different HPV types gene

Strain- 18	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus isolate HPV18	E1	KY501976.1	100	100	100
Human papillomavirus isolate HPV18	E1	KY501973.1	100	100	100
Human papillomavirus isolate HPV18 P1-20	E1	KY501967.1	100	100	100
Human papillomavirus isolate HPV18 P1	E1	KY501965.1	100	100	100
Human papillomavirus isolate HPV18 P2	E1	KY501971.1	100	100	100
Human papillomavirus isolate HPV18 P1-50	E1	KY501970.1	100	100	100
Human papillomavirus isolate HPV18 P2-30	E1	KY501974.1	100	100	100
Human papillomavirus isolate HPV18 P1-40	E1	KY501969.1	100	100	100

Stark Genotyping HPV Molecular Diagnostic Kit

CARBONTECHNOLOGIES

This analysis confirms the comprehensive coverage of the primer/probe sets across a diverse range of HPV sequences, ensuring the robust analytical sensitivity of the diagnostic Kit.

Strain- 45	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 45 isolate LNS8825361_HP45	E1	LR862061.1.1	100	100	100
Human papillomavirus type 45 isolate Qv34163	E1	KC470260.1	100	100	100
Human papillomavirus type 45 isolate Qv30712	E1	>KC470259.1.1	100	100	100
Human papillomavirus type 45 isolate Qv20214	E1	EF202156.1	100	100	100
Human papillomavirus type 45 isolate LNS8825361_HP45	E1	LR862061.1.1	100	100	100
Human papillomavirus type 45 isolate Qv27565	E1	EF202157.1	100	100	100
Human papillomavirus type 45 isolate Qv33330	E1	EF202158.1	100	100	100
Human papillomavirus type 45 isolate Qv34178	E1	EF202159.1	100	100	100

Strain- 31	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus 31 isolate VBD13/14	E2	MW814876.1	100	100	100
Human papillomavirus type 31 isolate LNS8465006_HP31	E2	LR862053.1	100	100	100
Human papillomavirus type 31 isolate LNS8357921_HP31	E2	LR862049.1	100	100	100
Human papillomavirus type 31 isolate LNS7548544_HP31	E2	LR862026.1	100	100	100
Human papillomavirus type 31 isolate LNS7384732	E2	LR862018.1	100	100	100
Human papillomavirus type 31 isolate LNS7358029	E2	LR862015.1	100	100	100
Human papillomavirus type 31 isolate LNS6001593	E2	LR861970.1	100	100	100
Human papillomavirus type 31 isolate LNS0937074	E2	LR861951.1	100	100	100
Human papillomavirus type 31 isolate IARC366181RW	E2	MT752571.1	100	100	100
Human papillomavirus type 31	E2	MT752570.1	100	100	100

Stark Genotyping HPV Molecular Diagnostic Kit

CARBONTECHNOLOGIES

isolate IARC366076RW					
Human papillomavirus type 31 isolate PAP293744	E2	MT752505.1	100	100	100
Human papillomavirus type 31 isolate PAP287325	E2	MT752449.1	100	100	100

Strain- 35	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 35 isolate LNS0838993	E6/E7	LR861946.1	100	100	100
Human papillomavirus type 35 isolate PAP294578	E6/E7	MT217995.1	100	100	100
Human papillomavirus type 35 isolate PAP282649	E6/E7	MT217962.1	100	100	100
Human papillomavirus type 35 isolate PAP163773	E6/E7	MT217679.1	100	100	100
Human papillomavirus type 35 isolate PAP162942	E6/E7	MT217678.1	100	100	100
Human papillomavirus type 35 isolate PAP145045	E6/E7	MT217620.1	100	100	100
Human papillomavirus type 35 isolate IARC1220350CH	E6/E7	MT217487.1	100	100	100
Human papillomavirus type 35 isolate IARC1101761NI	E6/E7	MT217470.1	100	100	100
Human papillomavirus type 35 isolate IARC1251067MO	E6/E7	MT217441.1	100	100	100
Human papillomavirus type 35 isolate IARC166070014RW	E6/E7	MT217402.1	100	100	100
Human papillomavirus type 35 isolate IARC1110712GU	E6/E7	MT217344.1	100	100	100
Human papillomavirus type 35 isolate IARC311027GU	E6/E7	MT217305.1	100	100	100

Strain- 59	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus 59 strain kyd-s0359	E6	MT783417.1	100	100	100
Human papillomavirus type 59 isolate LNS9786324	E6	LR862080.1	100	100	100
Human papillomavirus type 59 isolate LNS3279856	E6	LR861868.1	100	100	100
Human papillomavirus type 59	E6	KC470263.1	100	100	100

Stark Genotyping HPV Molecular Diagnostic Kit

CARBONTECHNOLOGIES

isolate Qv23880					
Human papillomavirus type 59 isolate LZod68-59	E6	EU918767.1	100	100	100
Human papillomavirus type 59 isolate LNS5859951	E6	LR861933.1	100	100	100
Human papilloma virus type 59	E6	X77858.1	100	100	100
Human papillomavirus type 59 isolate Qv25652	E6	KC470261.1	100	100	100
Human papillomavirus type 59 isolate Qv25808	E6	KC470264.1	100	100	100
Human papillomavirus type 59 isolate Qv33993	E6	KC470265.1	100	100	100
Human papillomavirus type 59 isolate 32A.59	E6	KU298921.1	100	100	100
Human papillomavirus type 59 isolate LNS2917513_HPV59	E6	LR861860.1	100	100	100

Strain- 39	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 39 isolate LNS9068892	L1	LR862071.1	100	100	100
Human papillomavirus type 39 isolate LNS8643648	L1	LR862054.1	100	100	100
Human papillomavirus type 39 isolate LNS7329845	L1	LR862014.1	100	100	100
Human papillomavirus type 39 isolate 39PL03 L1	L1	MK344660.1	100	100	100
Human papillomavirus type 39 isolate 39PL02 L1	L1	MK344659.1	100	100	100
Human papillomavirus type 39 isolate 16B	L1	KX514417.1	100	100	100
Human papillomavirus type 39 isolate Tw562	L1	KC470243.1	100	100	100
Human papillomavirus type 39 isolate Qv25609	L1	KC470242.1	100	100	100
Human papillomavirus type 39 isolate Qv25959	L1	KC470240.1	100	100	100
Human papillomavirus type 39 isolate Qv27715	L1	KC470237.1	100	100	100
Human papillomavirus type 39 isolate LNS4492831	L1	LR861904.1	100	100	100
Human papillomavirus type 39 clone 39L1.A major capsid	L1	JN104068.1	100	100	100

protein					
---------	--	--	--	--	--

Strain- 6/11	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 11 isolate HPV11-gw-1108	L1	MK463919.1	100	100	100
Human papillomavirus type 11 isolate HPV11-gw-1111	L1	MK463922.1	100	100	100
Human papillomavirus type 11 isolate HPV11-gw-1110	L1	MK463921.1	100	100	100
Human papillomavirus type 11 isolate HPV11-gw-1109	L1	MK463920.1	100	100	100
Human papillomavirus type 6 isolate HPV6-gw-0611	L1	MK463909.1	100	100	100
Human papillomavirus type 6 isolate HPV6-gw-0602	L1	MK463905.1	100	100	100
Human papillomavirus type 6 isolate ASCUS1/HPV6	L1	MK313781.1	100	100	100
Human papillomavirus type 6 isolate CAC1/HPV6	L1	MK313778.1	100	100	100
Human papillomavirus type 6 isolate JO-RRP2/HPV6	L1	MK313777.1	100	100	100
Human papillomavirus type 6 isolate AO-RRP7/HPV6	L1	MK313775.1	100	100	100
Human papillomavirus type 6 isolate AO-RRP5/HPV6	L1	MK313773.1	100	100	100
Human papillomavirus type 6 isolate AO-RRP3/HPV6	L1	MK313771.1	100	100	100
Human papillomavirus type 6 isolate AO-RRP1/HPV6	L1	MK313769.1	100	100	100
Human papillomavirus type 6 isolate 111B	L1	KX514423.1	100	100	100
Human papillomavirus type 6 isolate J-50	L1	KU721785.1	100	100	100

Strain- 51	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 51 isolate LNS9024319	L1	LR862069.1	100	100	100
Human papillomavirus type 51 isolate LNS8077695	L1	LR862046.1	100	100	100

Stark Genotyping HPV Molecular Diagnostic Kit

CARBONTECHNOLOGIES

Human papillomavirus type 51 isolate LNS0570411_HP51	L1	LR861939.1	100	100	100
Human papillomavirus type 51 isolate LNS4020261	L1	LR861893.1	100	100	100
Human papillomavirus type 51 isolate LNS3425544	L1	LR861869.1	100	100	100
Human papillomavirus type 51 isolate LNS9888540_HP51	L1	LR861810.1	100	100	100
Human papillomavirus type 51 isolate CL17 L1	L1	MH577961.1	100	100	100
Human papillomavirus type 51 isolate CL16	L1	MH577960.1	100	100	100
Human papillomavirus type 51 isolate CL10	L1	MH577954.1	100	100	100
Human papillomavirus type 51 isolate 118A.51	L1	KU298905.1	100	100	100

Strain- 66	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 66 isolate NGSk294-66	E6	LR861964.1	100	100	100
Human papillomavirus type 66 isolate LNS4474748_HP66	E6	LR861902.1	100	100	100
Human papillomavirus type 66 isolate LNS2720744	E6	LR861851.1	100	100	100
Human papillomavirus type 66 isolate 110	E6	JN661435.1	100	100	100
Human papillomavirus type 66 isolate 100	E6	JN661426.1	100	100	100
Human papillomavirus type 66 isolate 95	E6	JN661421.1	100	100	100
Human papillomavirus type 66 isolate 92	E6	JN661418.1	100	100	100
Human papillomavirus type 66 isolate 90	E6	JN661416.1	100	100	100
Human papillomavirus type 66 isolate Bsb-107	E6	HM585479.1	100	100	100
Human papillomavirus type 66 clone Qv11088	E6	EF177190.1	100	100	100

Strain- 68	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
------------	--------	-----------	---------------------------------	---------------------------------	------------------------

Human papillomavirus type 68 isolate 55-86-21	E6/E7	MK874671.1	100	100	100
Human papillomavirus type 68 isolate 21-54-34	E6/E7	MK874669.1	100	100	100
Human papillomavirus type 68 isolate 21-54-24	E6/E7	MK874668.1	100	100	100
Human papillomavirus type 68 isolate 21-54-18	E6/E7	MK874666.1	100	100	100
Human papillomavirus type 68 isolate 1-20-9	E6/E7	MK874662.1	100	100	100
Human papillomavirus type 68 isolate M180	E6/E7	MN047810.1	100	100	100
Human papillomavirus type 68 isolate 268	E6/E7	KX514428.1	100	100	100
Human papillomavirus type 68 isolate P141607	E6/E7	KU195243.1	100	100	100
Human papillomavirus type 68 isolate Rw826	E6/E7	KC470281.1	100	100	100
Human papillomavirus type 68 isolate Qv25395	E6/E7	KC470279.1	100	100	100

Strain- 56/66	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 56 isolate 56SE-23	E6/E7	KX645764.1	100	100	100
Human papillomavirus type 56 isolate 56SE-21	E6/E7	KX645762.1	100	100	100
Human papillomavirus type 56 isolate 56SE-19	E6/E7	KX645760.1	100	100	100
Human papillomavirus type 56 isolate 56SE-09	E6/E7	KX645750.1	100	100	100
Human papillomavirus type 56 isolate 60B	E6/E7	KX514418.1	100	100	100
Human papillomavirus type 56 isolate 87A.56	E6/E7	KU298916.1	100	100	100
Human papillomavirus type 56 clone Qv22608	E6/E7	EF177179.1	100	100	100
Human papillomavirus type 66 isolate LNS7384732	E6/E7	LR862020.1	100	100	100
Human papillomavirus type 66 isolate LNS1839688	E6/E7	LR861964.1	100	100	100
Human papillomavirus type 66 isolate LNS5491243	E6/E7	LR861922.1	100	100	100

Human papillomavirus type 66 isolate LNS3785675	E6/E7	LR861879.1	100	100	100
Human papillomavirus type 66 isolate 106	E6/E7	JN661479.1	100	100	100
Human papillomavirus type 66 isolate 95	E6/E7	JN661470.1	100	100	100
Human papillomavirus type 66 isolate 85	E6/E7	JN661460.1	100	100	100
Human papillomavirus type 66 isolate Qv32783	E6/E7	EF546480.1	100	100	100
Human papillomavirus type 66 clone Qv11088	E6/E7	EF177190.1	100	100	100

Strain- 33/52/ 58	Target	Accession	% Homology Test Forward primer%	% Homology Test Reverse primer%	% Homology Test Probe%
Human papillomavirus type 33 isolate LNS9453833	E1/E2	LR862077.1	100	100	100
Human papillomavirus type 33 isolate LNS3510974	E1/E2	LR861872.1	100	100	100
Human papillomavirus type 33 isolate 67B.33	E1/E2	KU298894.1	100	100	100
Human papillomavirus type 33 isolate 65C.33	E1/E2	KU298893.1	100	100	100
Human papillomavirus type 33 isolate BF375	E1/E2	KF436865.1	100	100	100
Human papillomavirus type 33	E1/E2	M12732.1	100	100	100
Human papillomavirus type 33 isolate LZcc12-33	E1/E2	EU918766.1	100	100	100
Human papillomavirus type 33 isolate Qv32494	E1/E2	HQ537688.1	100	100	100
Human papillomavirus 52 isolate PAP3497	E1/E2	MZ374448.1	100	100	100
Human papillomavirus 52 isolate PAP1639	E1/E2	MZ374436.1	100	100	100
Human papillomavirus 52 isolate PAP1189	E1/E2	MZ374435.1	100	100	100
Human papillomavirus 52 isolate PAP2912	E1/E2	MZ374431.1	100	100	100
Human papillomavirus 52 isolate PAP2357	E1/E2	MZ374425.1	100	100	100
Human papillomavirus 52 isolate PAP2484	E1/E2	MZ374419.1	100	100	100
Human papillomavirus 52	E1/E2	MZ374415.1	100	100	100

isolate PAP1917					
Human papillomavirus 52 isolate PAP2698	E1/E2	MZ374413.1	100	100	100
Human papillomavirus type 52 isolate LNS9453833	E1/E2	LR862077.1	100	100	100
Human papillomavirus type 52 isolate LNS6759684	E1/E2	LR861998.1	100	100	100
Human papillomavirus type 58 isolate LNS6048641	E1/E2	LR861973.1	100	100	100
Human papillomavirus type 58 isolate LNS0838993	E1/E2	LR861950.1	100	100	100
Human papillomavirus type 58 isolate LNS5390598	E1/E2	LR861919.1	100	100	100
Human papillomavirus type 58 isolate LNS4582218	E1/E2	LR861906.1	100	100	100

As supplemental data, wet testing of inclusivity was conducted using the DNall VirAll Kit. Three HPV-positive specimens from the WHO HPV Lab Net were tested, all of which were confirmed positive by the Stark Genotyping HPV Molecular Diagnostic Kit. Each specimen was diluted to $\leq 3\log_{10}$ LOD, $\leq 2\log_{10}$ LOD, and $\leq 1\log_{10}$ LOD in a negative specimen matrix (Cervical swab) and tested in the tenth replicate. The table below provides details of this testing:

This evaluation confirms the clinical sensitivity of the Stark Genotyping HPV Molecular Diagnostic Kit across a range of dilutions in a negative specimen matrix, demonstrating its ability to accurately detect HPV-positive specimens.

Table 9: the results of diagnostic sensitivity of Stark Quantitative HPV Molecular Diagnostic Kit by CARBON Technologies LLC

Dilution series	Copies/5 μ l	Ct
3log10 > LOD	50.000	24.55
		24.25
		25.05
		25.13
		24.71
		24.87
		24.18
		24.46
		25.14
		25.23
2log10 > LOD	5.000	28.44
		28.48
		29.31
		28.64
		28.78
		28.73
		29.08

		29.17
		29.08
		28.87
1log10 > LOD	500	34.84
		33.82
		34.52
		33.88
		34.17
		35.27
		35.37
		35.76
		35.78
		35.37

- Cross-reactivity (Analytical Specificity)**

The cross-reactivity of the Stark Genotyping HPV Molecular Diagnostic Kit was evaluated through in silico analysis and wet testing of potentially cross-reactive whole pathogens or purified nucleic acid from clinical specimens. No cross-reactivity was detected during either analysis.

In-silico mapping analysis of each primer/probe against several pathogens was conducted using the NCBI nr/nt database, accessed on January 9, 2022, through the online BLASTN 2.10.0+ tool. Representative results of this analysis are provided below the table.

No cross-reactivity was observed for the listed pathogens in either the in-silico or wet-testing experiments, confirming the high analytical specificity of the diagnostic Kit.

Table 10: The In-Silico Specificity Analysis of Primer and Probe Set for Other pathogens.

Strain- 18 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	HIV-1 isolate 002211_C07	MT033226.1	59%	72%	62%
Hepatitis C virus subtype 1a	10jsszGP017	JQ303503.1	63%	55%	62%
HSV-1	s17pp22a	MN159382.1	70%	50%	66%
Hepatitis B virus	isolate 2793	KF169309.1	50%	66%	59%
HSV-2	strain 2009-24855_S55_L001	MH790632.1	45%	61%	70%
Mycoplasma genitalium	M2321 MgPar 9	CP003773.1	45%	55%	51%
Chlamydia trachomatis	TC0350	CP042803.1	66%	61%	77%
Streptococcus agalactiae	Sag27	CP031556.1	81%	77%	48%
Human T-cell leukemia virus	057-OL	MN453165.1	50%	88%	68%

type I					
Human gamma herpesvirus	HC-0036	LR813015.1	63%	63%	62%
Human alpha herpesvirus	s17pp22a	MN159382.1	95%	55%	66%
Human T-lymphotropic virus 2	HTLV-2/Japan/NIID/001/2020	K01670.1	50%	44%	51%
Human parvovirus B19	911NewF	MN088124.1	33%	55%	40%
JC polyomavirus	UJC005	MW587996.1	72%	50%	88%
Neisseria gonorrhoeae	TUM16691	AP023075.1	54%	55%	44%
Trichomonas vaginalis	TVAG_290210	XM_001312003.1	59%	72%	62%

Strain- 16 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	YBF319	KT198590.1	52%	63%	65%
Hepatitis C virus subtype 1a	HCV-1b/US/BID-V377/2006	EU256065.1	48%	40%	55%
HSV-1	2007-16123	MG999881.1	36%	33%	50%
Hepatitis B virus	B127W	MH051987.1	48%	40%	50%
HSV-2	2006-15095	MH899846.1	55%	40%	72%
Mycoplasma genitalium	1206KLM26 MgpC adhesin	MT439445.1	44%	53%	50%
Chlamydia trachomatis	TC0411:TncL2/tetR-104	CP042798.1	75%	86%	75%
Streptococcus agalactiae	2013-1366	CP051844.1	70%	84%	44%
Human T-cell leukemia virus type I	Aus-CS	KF242506.1	60%	45%	30%
Human gammaherpesvirus	21_AJOUA	MT648662.1	55%	48%	43%
Human alphaherpesvirus	HSV1-ORIGINAL-H1	MH160367.1	36%	33%	50%
Human T-lymphotropic virus 2	(HTLV-II AA	X76030.1	40%	72%	43%

Human parvovirus B19	HZ13 NS1	KT389459.1	43%	73%	75%
JC polyomavirus	IRN-26H	MN163027.1	70%	36%	66%
Neisseria gonorrhoeae	NG196	CP043872.1	40%	40%	44%
Trichomonas vaginalis	TVAG_231320	XM_001313720.1	52%	46%	60%

Strain- 31 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	A09.Q2.20M9D1	MH264318.1	63%	63%	68%
Hepatitis C virus subtype 1a	AV67	EU484132.1	58%	75%	52%
HSV-1	HSV-H1412	MH999851.1	56%	87%	45%
Hepatitis B virus	551.31	KR013975.1	70%	58%	52%
HSV-2	HSV2_19	MH697440.1	64%	87%	52%
Mycoplasma genitalium	1493B26 MgpB	MT439517.1	52%	56%	42%
Chlamydia trachomatis	TC0657:TnxL2/tetR-403	CP042757.1	87%	42%	54%
Streptococcus agalactiae	CM45/1	CP088933.1	55%	87%	62%
Human T-cell leukemia virus type I	IR (94) bZIP	MN453136.1	57%	47%	52%
Human gammaherpesvirus	20_IHOM	MT648661.1	78%	48%	58%
Human alphaherpesvirus	CM1v8	KX791821.1	56%	87%	56%
Human T-lymphotropic virus 2	N1232	KY493238.1	50%	75%	42%
Human parvovirus B19	NA121	MG765334.1	50%	56%	78%
JC polyomavirus	SA38527_04 VP1	EU835186.1	42%	42%	87%
Neisseria gonorrhoeae	SRRSH205	CP048907.1	78%	68%	62%
Trichomonas vaginalis	TVAG_449190	XM_001308315.1	62%	63%	68%

Strain- 35	Strain	GenBank Acc#	%	%	%
------------	--------	--------------	---	---	---

Stark Genotyping HPV Molecular Diagnostic Kit

CARBONTECHNOLOGIES

Pathogen (Taxonomy ID)			Homology Test FP	Homology Test RP	Homology Test Probe
Human immunodeficiency virus	OU072	OK180789.1	48%	44%	56%
Hepatitis C virus subtype 1a	EA03T0	AM271169.1	47%	39%	43%
HSV-1	2011-12741	MG999893.1	43%	74%	47%
Hepatitis B virus	CHB-47-PS	KY428722.1	40%	37%	43%
HSV-2	2012-5447	MH790637.1	56%	43%	52%
Mycoplasma genitalium	1491KLM26	MT439589.1	33%	52%	44%
Chlamydia trachomatis	TC0411: TnxD/tetR-103	CP042802.1	54%	81%	48%
Streptococcus agalactiae	SA627	CP019837.1	43%	52%	33%
Human T-cell leukemia virus type I	BHP00140	MG388047.1	52%	43%	44%
Human gammaherpesvirus	19_ELTW	MT648660.1	40%	48%	47%
Human alphaherpesvirus	2011-17239	MG999896.1	43%	74%	47%
Human T-lymphotropic virus 2	BBD_7003	FJ911659.1	51%	29%	59%
Human parvovirus B19	D1599 NS1	DQ234775.2	40%	55%	59%
JC polyomavirus	J007 VP1	MK477564.1	37%	88%	33%
Neisseria gonorrhoeae	TUM19853	AP023067.1	47%	77%	52%
Trichomonas vaginalis	TVAG_196200	XM_001312988.1	51%	66%	52%

Strain- 33/52/58 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	AA075a_WG7	JX447542.1	63%	66%	59%

Hepatitis C virus subtype 1a	PATNO72	AB285080.1	45%	62%	41%
HSV-1	2004-63623	MG999877.1	62%	59%	72%
Hepatitis B virus	OHBV-HIV137	MF618348.1	45%	45%	50%
HSV-2	2006-18003CAM	MH790593.1	41%	59%	62%
Mycoplasma genitalium	10467 MgpB	KP318824.1	50%	63%	40%
Chlamydia trachomatis	TC0189:TncL2/tetR-102	CP042734.1	87%	59%	63%
Streptococcus agalactiae	NEM316	AL766847.1	45%	87%	50%
Human T-cell leukemia virus type I	BHP00005	MG388044.1	75%	59%	54%
Human gammaherpesvirus	sLCL-T3.27	LC573551.1	59%	50%	45%
Human alphaherpesvirus	2003-15756	MG999872.1	75%	68%	62%
Human T-lymphotropic virus 2	NIID18001	LC440555.1	36%	41%	31%
Human parvovirus B19	BN58.3	DQ408302.1	54%	77%	66%
JC polyomavirus	IRN-72H	MN163031.1	45%	37%	63%
Neisseria gonorrhoeae	SS3160	AP019853.2	41%	58%	45%
Trichomonas vaginalis	TVAG_308110	XM_001320393.1	40%	58%	72%

Strain- 59 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	505_1982a	MG197155.1	60%	52%	68%
Hepatitis C virus subtype 1a	GZ52540	KC844051.1	35%	42%	46%
HSV-1	11_DOCK8	MN401208.1	64%	40%	42%
Hepatitis B virus	NOA_142	MG098582.1	48%	40%	44%
HSV-2	HJ12	MN187895.1	40%	36%	52%
Mycoplasma genitalium	M2288	CP003773.1	44%	56%	40%
Chlamydia	TC0350:TncL2/tetR-936	CP042803.1	44%	40%	64%

trachomatis					
Streptococcus agalactiae	NEM316	AL766845.1	56%	40%	36%
Human T-cell leukemia virus type I	PH1511	MN781156.1	40%	56%	40%
Human gammaherpesvirus	HKHD138	MH590507.1	64%	44%	76%
Human alphaherpesvirus	2007-16123	MG999881.1	64%	40%	36%
Human T-lymphotropic virus 2	BBD_3079	FJ911654.1	36%	40%	32%
Human parvovirus B19	HZ13 NS1	KT389459.1	40%	36%	76%
JC polyomavirus	Tn-19 VP1	JQ433658.1	36%	32%	48%
Neisseria gonorrhoeae	TUM16691	AP023075.1	44%	76%	52%
Trichomonas vaginalis	TVAG_149710	XM_001308182.1	60%	52%	56%

Strain- 39 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	505_1982a	MG197155.1	60%	52%	68%
Hepatitis C virus subtype 1a	GZ52540	KC844051.1	35%	42%	46%
HSV-1	11_DOCK8	MN401208.1	64%	40%	42%
Hepatitis B virus	NOA_142	MG098582.1	48%	40%	44%
HSV-2	HJ12	MN187895.1	40%	36%	52%
Mycoplasma genitalium	M2288	CP003773.1	44%	56%	40%
Chlamydia trachomatis	TC0350:TncL2/tetR-936	CP042803.1	44%	40%	64%
Streptococcus agalactiae	NEM316	AL766845.1	56%	40%	36%
Human T-cell leukemia virus type I	PH1511	MN781156.1	40%	56%	40%
Human gammaherpesvirus	HKHD138	MH590507.1	64%	44%	76%
Human alphaherpesvirus	2007-16123	MG999881.1	64%	40%	36%
Human	BBD_3079	FJ911654.1	36%	40%	32%

T-lymphotropic virus 2					
Human parvovirus B19	HZ13 NS1	KT389459.1	40%	36%	76%
JC polyomavirus	Tn-19 VP1	JQ433658.1	36%	32%	48%
Neisseria gonorrhoeae	TUM16691	AP023075.1	44%	76%	52%
Trichomonas vaginalis	TVAG_149710	XM_001308182.1	60%	52%	56%

Strain- 66 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	2303-08PT46	MH832382.1	50%	35%	46%
Hepatitis C virus subtype 1a	EG06T1ak	AM708703.1	35%	42%	39%
HSV-1	H166	KM222726.1	35%	39%	53%
Hepatitis B virus	124_CA_Kal	MH580631.1	39%	35%	42%
HSV-2	HJ12	MN187895.1	35%	50%	57%
Mycoplasma genitalium	199 MgPar 7	EF117299.2	57%	42%	39%
Chlamydia trachomatis	KU2043c768-137	CP042789.1	32%	39%	35%
Streptococcus agalactiae	NY84115	JF270524.1	35%	39%	42%
Human T-cell leukemia virus type I	Gab1014FC	EU444097.1	80%	50%	38%
Human gammaherpesvirus	DOCK8	MN401207.1	35%	42%	53%
Human alphaherpesvirus	E25	HM585506.2	50%	57%	42%
Human T-lymphotropic virus 2	HTLV-2/Japan/NIID/001/2020	LC534557.1	57%	50%	73%
Human parvovirus B19	BX2	MH201456.1	34%	30%	65%
JC polyomavirus	VP1	EF369493.1	61%	69%	42%
Neisseria gonorrhoeae	SW0236	CP061485.1	61%	53%	30%
Trichomonas vaginalis	TVAG_402520	XM_001580933.1	42%	69%	50%

Strain- 68 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	BA855	MK732623.1	57%	76%	50%
Hepatitis C virus subtype 1a	9419KBS/12	MH627077.1	50%	53%	50%
HSV-1	DOCK8	MN401208.1	46%	61%	73%
Hepatitis B virus	BB1618958	MH247244.1	69%	42%	57%
HSV-2	HJ12	MN187895.1	61%		
Mycoplasma genitalium	TC0350:TncL2/tetR-936	CP042803.1	35%	39%	57%
Chlamydia trachomatis	01173	CP053027.1	42%	69%	42%
Streptococcus agalactiae	IL1657	KF202322.1	61%	36%	34%
Human T-cell leukemia virus type I	IT1036	KF202320.1	34%	32%	38%
Human gammaherpesvirus	HKHD44	MH590413.1	42%	38%	69%
Human alphaherpesvirus	BB1618958_p	MH247244.1	30%	61%	50%
Human T-lymphotropic virus 2	HTLV-2/Japan/NIID/001/2020	LC534557.1	69%	34%	30%
Human parvovirus B19	R0416	DQ234769.2	38%	50%	42%
JC polyomavirus	JCV144CSFRRC	JF425684.1	65%	46%	57%
Neisseria gonorrhoeae	TUM16691	AP023075.1	42%	65%	46%
Trichomonas vaginalis	TVAG_454270	XM_001329699.1	46%	50%	42%

Strain- 6/11 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	7864_6	MW189577.1	54%	41%	50%
Hepatitis C virus subtype 1a	1b.NR.022	KM580649.1	66%	50%	54%

HSV-1	Sc16	MN159383.1	54%	45%	50%
Hepatitis B virus	cxx16	MZ244221.1	50%	41%	
HSV-2	HJ12	MN187895.1	45%	45%	54%
Mycoplasma genitalium	M2288	CP003773.1	62%	54%	41%
Chlamydia trachomatis	TC0693:TncL2/tetR-3509	CP042693.1	50%		41%
Streptococcus agalactiae	QMA0271	CP029632.1	50%	54%	66%
Human T-cell leukemia virus type I	mel5	L02534.1	54%	41%	45%
Human gammaherpesvirus	RK_LCL_L3	MG298914.1	66%	45%	54%
Human alphaherpesvirus	McKrae	MN136524.1	37%	50%	70%
Human T-lymphotropic virus 2	HTLV-2/Japan/NIID/001/2020	LC534557.1	70%	37%	66%
Human parvovirus B19	057214	JN211130.1	29%	62%	45%
JC polyomavirus	NIID11-53	LC164352.1	70%	33%	37%
Neisseria gonorrhoeae	TUM16691	AP023075.1	66%	37%	45%
Trichomonas vaginalis	TVAG_439770	XM_001306803.1	62%	29%	37%

Strain- 56/66 Pathogen (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	RF (HAT-3)	M17451.1	83%	50%	58%
Hepatitis C virus subtype 1a	Pt9.Seg8-Con	KX084702.1	50%	75%	34%
HSV-1	2339/2009	JQ352258.1	45%	70%	66%
Hepatitis B virus	HBV_NWFD_AHB111	MZ090870.1	75%	54%	58%
HSV-2	M2288	CP003773.1	58%	50%	
Mycoplasma genitalium	TC0411: TnxD/tetR-103	CP042802.1	75%	62%	58%
Chlamydia trachomatis	TC0411: TnxD/tetR-103	CP042802.1	50%	75%	45%

Streptococcus agalactiae	SG-M25	CP021867.1	58%	54%	41%
Human T-cell leukemia virus type I	AC038	LC192501.1	41%	54%	58%
Human gammaherpesvirus	eBL-Tumor-0014	LR813032.1	30%	41%	
Human alphaherpesvirus	2339/2009	JQ352258.1	62%	56%	50%
Human T-lymphotropic virus 2	HTLV-2/Japan/NIID/001/2020	LC534557.1	41%	50%	45%
Human parvovirus B19	HZ92	KT389469.1	34%	30%	58%
JC polyomavirus	RKAB15 LTag	LC650363.1	35%	42%	61%
Neisseria gonorrhoeae	FA6140	CP012027.1	34%	62%	58%
Trichomonas vaginalis	TVAG_480160	XM_001307664.1	50%	30%	62%

Pathogen / 51 (Taxonomy ID)	Strain	GenBank Acc#	% Homology Test FP	% Homology Test RP	% Homology Test Probe
Human immunodeficiency virus	20263v08	MN791407.1	54%	65%	50%
Hepatitis C virus subtype 1a	HCV-1b/US/BID-V152/2003	EU155224.2	45%	50%	50%
Hepatitis B virus	94T.13.	KR812069.1	54%	50%	54%
HSV-1	2018-5971	MH813987.1	87%	45%	59%
Human papillomavirus	NCI_230248	MG849960.1	40%	45%	45%
HSV-2	1996-45091	MH790670.1	66%	83%	45%
Mycoplasma genitalium	M6285 MgPa	GU226202.1	59%	45%	41%
Chlamydia trachomatis	TC0411: TnxD/tetR-296	CP042800.1	50%	63%	81%
Streptococcus agalactiae	SG-M25	CP021867.1	68%	65%	83%
Human T-cell leukemia virus type I	012BR_HAM107	KF797884.1	41%	50%	40%
Human gammaherpesvirus	HC-0020	LR813006.1	54%	85%	58%

Human alphaherpesvirus	Y033 UL40	AY240581.1	45%	87%	50%
Human T-lymphotropic virus 2	NIID18002	LC440556.1	40%	45%	45%
Human parvovirus B19	CAEN.FRA 21.09/1	FN669504.1	37%	40%	45%
JC polyomavirus	PE-21	AB081027.1	33%	45%	40%
Neisseria gonorrhoeae	AUSMDU00010541	CP045832.1	63%	55%	45%
Trichomonas vaginalis	TVAG_308030	XM_001579636.1	59%	55%	50%

- **Cross-Reactivity (Clinical Specificity)** of the Stark Genotyping HPV Molecular Diagnostic Kit

To assess the clinical specificity of the Kit, nucleic acid from other pathogens was tested in a matrix of negative samples (Cervical negative swab) diluted to a specific concentration. These samples were then extracted and tested using the Stark Genotyping HPV Molecular Diagnostic Kit. No cross-reactivity was observed for the pathogens listed in the following table:

This testing confirms the high clinical specificity of the Stark Genotyping HPV Molecular Diagnostic Kit, indicating its ability to accurately detect HPV without interference from other pathogens.

Table 12: Investigation of the cross-reactivity of the HPV using Stark Genotyping HPV Molecular Diagnostic Kit

Virus/Bacteria/Parasite	Source/ Sample type	Ct Value
Human immunodeficiency virus-1	Clinical sample	-/-
Hepatitis C virus	Clinical sample	-/-
Cytomegalovirus	Clinical sample	-/-
Herpes simplex virus type 1	Clinical sample	-/-
Herpes simplex virus type 2	Clinical sample	-/-
Human papillomavirus	Clinical sample	-/-
Epstein-Barr virus	Clinical sample	-/-
Adenovirus	Clinical sample	-/-
Influenza A	Clinical sample	-/-
Influenza B	Clinical sample	-/-
Legionella pneumophila	Clinical sample	-/-
Cryptococcus neoformans	Clinical sample	-/-
Chlamydia pneumonia	Clinical sample	-/-
Streptococcus pneumonia	Clinical sample	-/-
Respiratory Syncytial Virus	Clinical sample	-/-
Mycoplasma pneumonia	Clinical sample	-/-
Streptococcus pyogenes	Clinical sample	-/-
Mycobacterium tuberculosis	Clinical sample	-/-
10 Pooled human genomes	Clinical sample	-/-

- **Clinical Evaluation**

Stark Genotyping HPV Molecular Diagnostic Kit

The clinical performance of the Stark Genotyping HPV Molecular Diagnostic Kit was assessed using 206 cervical specimens, liquid-based cytology specimens, urine, and paraffin-embedded tissue collected from patients suspected of HPV infection. The comparator methods utilized were the HPV 3.5 LCD-Array Kits (Chipron GmbH) and the Cobas® HPV test (Roche), both of which have received CE-IVD certification. The extraction method employed was the DNall VirAll Kit.

The results of the clinical evaluation are summarized in the analysis, demonstrating a Positive Percent Agreement (PPA) of 96.47% and Negative Percent Agreement (NPA) of 100% when compared to the HPV 3.5 LCD-Array Kits (Chipron GmbH). Similarly, when compared to the Cobas® HPV test (Roche), the Stark Genotyping HPV Molecular Diagnostic Kit exhibited a PPA of 94.4% and NPA of 100%.

Table 13: Results of the Stark Genotyping HPV Molecular Diagnostic Kit (CARBON Technologies LLC) alongside the HPV 3.5 LCD-Array Kits (Chipron GmbH) and Cobas® HPV test (Roche).

Number	HPV 3.5 LCD-Array Kits (Chipron GmbH)	Stark Genotyping HPV Molecular Diagnostic Kit (CARBON Technologies LLC)
	Result	Result
1	11	6/11
2	6/16	6/11, 16
3	51	-
4	16	16
5	51	51
6	6	6/11 33/52/58
7	6/59	6/11
8	16/18	16/18
9	16	16
10	45	45
11	31	31
12	45	45
13	45	45
14	45	45
15	18	18
16	18	18
17	18	18
18	45	45
19	45	45
20	16	16
21	16	16
22	16	16
23	16	16
24	16	16
25	16	16
26	16	16
27	51	51
28	51	51
29	51	51

30	51	51
31	51	51
32	56	56/66
33	56	56/66
34	56	56/66
35	56	56/66
36	56	56/66
37	56	56/66
38	56	56/66
39	66	56/66
40	66	56/66
41	66	56/66
42	66	56/66
43	66	56/66
44	66	56/66
45	66	56/66
46	35	35/39
47	35	35/39
48	35	35/39
49	35	35/39
50	35	35/39
51	35	35/39
52	35	35/39
53	35	-
54	59	59
55	59	59
56	59	59
57	59	59
58	59	59
59	59	59
60	59	-
61	6	6/11
62	6	6/11
63	11	6/11
64	11	6/11
65	11	6/11
66	52	33/52/58
67	52	33/52/58
68	52	33/52/58
69	52	33/52/58
70	52	33/52/58
71	52	33/52/58
72	58	33/52/58

73	58	33/52/58
74	58	33/52/58
75	58	33/52/58
76	58	33/52/58
77	58/33	33/52/58
78	58/33	33/52/58
79	33	33/52/58
80	33	33/52/58
81	33	33/52/58
82	33	33/52/58
Number	Cobas® HPV test (Roche)	Stark Genotyping HPV Molecular Diagnostic Kit (CARBON Technologies LLC)
1	16	16
2	16	16
3	16	16
4	16	16
5	16	16
6	16	16
7	16	16
8	16	16
9	16	16
10	16	16
11	16	16
12	16	-
13	18	18
14	18	18
15	18	18
16	18	18
17	18	18
18	18	18
19	18	18
20	18	18
21	18	18
22	18	-
23	High Risk	33/52/58, 56/66
24	High Risk	6/11, 33/52/58, 51
25	High Risk	6/11, 16, 51
26	High Risk	6/11, 51
27	High Risk	6/11, 33/52/58, 68, 51
28	High Risk	33/52/58, 51
29	High Risk	59
30	High Risk	33/52/58, 51

31	High Risk	33/52/58, 68, 51
32	High Risk	39, 33/52/58, 35, 68
33	High Risk	68, 51, 56/66
34	High Risk	51, 56/66

- **The results of Clinical Evaluation**

Table 14: Clinical Evaluation of Stark Genotyping HPV Molecular Diagnostic Kit compared to the HPV 3.5 LCD-Array Kits (Chipron GmbH)

Test		HPV 3.5 LCD-Array Kits (Chipron GmbH)		Total
		Positive	Negative	
Stark Genotyping HPV Molecular Diagnostic Kit (CARBON Technologies LLC)	Positive	82	0	82
	Negative	3	100	103
Total		85	100	185

- Positive Agreement Rate: $82 \div 85 \times 100\% = 96.47\%$
- Negative Agreement Rate: $100 \div 100 \times 100\% = 100\%$
- Overall rates of agreement: $(100 + 82) \div (82 + 0 + 100 + 3) \times 100\% = 98.37\%$

These rates indicate the level of agreement between the Stark Genotyping HPV Molecular Diagnostic Kit and the HPV 3.5 LCD-Array Kits by Chipron GmbH in detecting HPV infections, demonstrating high levels of concordance between the two diagnostic methods.













Table 15: Clinical Evaluation of Stark Genotyping HPV Molecular Diagnostic Kit compared to the Cobas® HPV test (Roche)

Test		Cobas® HPV test (Roche)		Total
		Positive	Negative	
Stark Genotyping HPV Molecular Diagnostic Kit (CARBON Technologies LLC)	Positive	34	0	34
	Negative	2	100	102
Total		36	100	136

- Positive Agreement Rate: $34 \div 36 \times 100\% = 94.44\%$
- Negative Agreement Rate: $100 \div 100 \times 100\% = 100\%$
- Overall rates of agreement: $(100 + 34) \div (34 + 0 + 100 + 2) \times 100\% = 98.52\%$

These rates indicate the level of agreement between the Stark Genotyping HPV Molecular Diagnostic Kit and the Cobas® HPV test by Roche in detecting HPV infections, demonstrating high levels of concordance between the two diagnostic methods.

Symbols

	Manufacturer
	Date of manufacture
	Use-by date
	Contains sufficient for <n> tests
	Consult instructions for use
	Biological risks
	Temperature limit
	CE Marking
	EU Representative
	In Vitro diagnostic medical device
	Catalogue Number
	Lot Number

Troubleshooting

For troubleshooting assistance with the Stark Genotyping HPV Molecular Diagnostic Kit, please refer to the following guidelines. Our ROJE Technical Support Team is available to address any further questions or concerns you may have:

Problems	Possible Causes	Action
No fluorescent signal is detected in any samples, including positive control.	Error in the preparation of the master mixture	Verify each component and ensure the volumes of reagent dispensed during the preparation of the master mixture are correct. Repeat PCR mixture preparation.
	Instrument settings error	Verify the Real-time PCR instrument settings are correct.
If the fluorescent signal is detected in a negative control reaction	Contamination of the extraction/preparation area	Clean surfaces and instruments with aqueous detergents, wash lab coats and replace test tubes and tips in use.
	PCR tube not properly sealed	Ensure plates are sealed correctly
If the fluorescent signal does not display the sigmoidal characteristic	Components degraded	Use a new batch.
	Poor quality of DNA samples carrying interferences	Repeat the test with the neat extracted DNA and 1:2 dilution of the extracted DNA.

	PCR equipment failure	Repeat the test or contact the equipment supplier
--	-----------------------	---

Technical assistance

For technical assistance, CARBON Technologies LLC guarantees your complete satisfaction. Our ROJE technical support team is comprised of highly trained and experienced scientists who are capable of troubleshooting most problems you may encounter. Our technical support team can offer expert advice to help you select the most suitable product for your needs.

You can contact our technical support team at any time through the following methods:

- By phone or fax: +96897058350
- Directly submit your questions to the CARBON Technologies technical support team through our website: www.carbontechnologiesco.com
- Send your questions via email to: technicalsupport@carbontechnologiesco.com

Rest assured, our team is dedicated to providing prompt and effective assistance to address any inquiries or issues you may have. We are committed to ensuring your satisfaction and success in using our products.

References

1-Ferlay J, Ervik M, Lam F, Colombet M, Mery L, Piñeros M, Znaor A, Soerjomataram I, Bray F (2018). Global Cancer Observatory: Cancer Today. Lyon, France: International Agency for Research on Cancer. Available from: <https://gco.iarc.fr/today>

2-Stelzle D, Tanaka LF, Lee KK, et al. Estimates of the global burden of cervical cancer associated with HIV. Lancet Glob Health 2020; published online Nov 16. DOI: S2214-109X (20)30459-9 [https://www.thelancet.com/journals/langlo/article/PIIS2214-109X\(20\)30459-9/fulltext](https://www.thelancet.com/journals/langlo/article/PIIS2214-109X(20)30459-9/fulltext)

3-Lei et al. (2020) HPV Vaccination and the Risk of Invasive Cervical Cancer. N Engl J Med 2020; 383:1340-8. DOI: 10.1056/NEJMoa1917338

4-World Health Organization. Global strategy to accelerate the elimination of cervical cancer as a public health problem <https://www.who.int/publications/i/item/9789240014107>

Contact Information



Medunion S.L.
Carrer de Tapioles 33, 2-1, 08004, Barcelona, SPAIN.



Carbon Technologies LLC
Innovation Park Muscat (IPM), P.O. Box 92, Al Khoudh 123, Muscat, OMAN.

24-hour service hotline: +968-97058350

After-sale Service Center: Carbon Technologies LLC



Release Date: Date of Manufacture: